

Image 1: Catheter Guide with NecroMax 6.0.

# EndoRotor<sup>®</sup> Catheter Guide<sup>™</sup>

The EndoRotor Catheter Guide device is intended to provide a 6mm auxiliary working channel for the EndoRotor Catheter during EndoRotor procedures.

## Power of Positioning

Catheter Guide can be positioned axially at different points on the scope during set-up. Endoscopists have the power to choose the catheters orientation during a procedure.




Image 2: EndoRotor positioned at 7 and 5 o'clock.

## Scope Compatibility


Catheter Guide is compatible with most diagnostic gastroscopes (reference IFU). The combined outer diameter when attached to a diagnostic gastroscope is 17-18mm.

Part Number	Description
6.0-PED-EGD	EndoRotor 6.0 PED Catheter Gastro Olympus
3.2-PED-EGD-OP	EndoRotor PED Catheter, Gastro Olympus
3.2-PED-EGD-F	EndoRotor PED Catheter, Gastro Fuji
CATHGUIDE	EndoRotor Catheter Guide
PED-PURGE	EndoRotor PED Purge Kit (Box 10)

For more information visit [www.interscopemed.com](http://www.interscopemed.com) or contact Customer Service:

### USA

Tel: +1-877-420-7299 (toll-free in U.S.)  
 Email: [info@interscopemed.com](mailto:info@interscopemed.com)

  
 Interscope, Inc.  
 200 Commerce Drive  
 Northbridge, MA 01534  
 USA  
 TEL: +1-800-461-4289  
 0031 (0) 73-3032541  
[www.interscopemed.com](http://www.interscopemed.com)

  
 EMERGO EUROPE  
 Prinsesgracht 20  
 2514 AP The Hague  
 The Netherlands  
[EMERGObyUL.com](http://EMERGObyUL.com)

  
 2797